

Mobile 12-Lead ECG

Quick Reference Guide

Kardia 12L is capable of recording a resting 12-lead ECG when connected to a compatible smart device (alivecor.com/compatibility/pro).

Use the QR code to download the KardiaStation app before use



Included components

- Kardia 12L ECG Module with (2) AAA Alkaline Batteries
- Patient Lead Wire
- Snap to Tab Adapters (5)
- Screwdriver for battery replacement

Visit alivecor.com/Kardia12L for complete Instructions for Use, and all Warnings and Precautions.


For additional information, or to request a paper copy of the Instructions for Use, contact support@alivecor.com or call +1 (855) 338-8800.




This set of lead wires meets the requirements of the standard, ECG trunk cables and patient lead wires (ANSI/AAMI EC53). This set of lead wires is intended for multi-patient use.

Lead placement options

Electrodes can be placed on the supine patient in either the Default or Alternate lead set locations. The KardiaStation app will guide you to select the appropriate lead placement, based on patient symptoms.



Default



Alternate

